



Solutions for better hearing

batteries

- superior quality
- longer lasting
- environmentally friendly





audifon Solutions for better hearing

audifon offers the complete range of hearing aid zinc air batteries. audifon hearing aid batteries have been tested and proven to be highly reliable and consistent in quality.

battery size	colour code	capacity
675	blue	640 mAh
13	orange	310 mAh
312	brown	180 mAh
10	yellow	105 mAh

Technological information

- Battery sizes available from 675 to 10.
- Battery power depends on amplification and usage period of the hearing system.
- Store at room temperature: heat and coldness may shorten the product life of the batteries.

Corporate Headquarters

audifon GmbH & Co. KG
 Werner-von-Siemens-Straße 2
 99625 Kölleda
 Germany

Phone +49-3635-4056-590
 Fax +49-3635-4056-589

contact@audifon.com

www.audifon.com