



Solutions for better hearing

looplink

- Bluetooth based accessory
- convenient hands-free communication for your mobile phone
- connect to all Bluetooth enabled devices easily





audifon Solutions for better hearing

The looplink is a lightweight wireless device utilising Bluetooth technology that provides hands-free communication, for example from a mobile phone to the hearing system.

This device is powered by AAA battery (alkaline or rechargeable). The looplink can be worn like a necklace.

Technological information

- looplink is compatible with devices using Bluetooth version 1.2/2.0
- signal radius area of 10 metres
- weight: 11,5 g (without battery)
- AAA batteries weigh depending on the producer 9 g–12 g
- power on standby mode: up to 500 hours
- talk time: up to 15 hours

Corporate Headquarters

audifon GmbH & Co. KG
Werner-von-Siemens-Straße 2
99625 Kölleda
Germany

Phone +49-3635-4056-590
Fax +49-3635-4056-589

contact@audifon.com

www.audifon.com