

Solutions for better hearing

**RadioLight DIR II**

- light weight and wireless headphones
- one transmitter for several headphones
- radio transmission up to 80 m





## audifon Solutions for better hearing

Ergonomical design and light weight, the wireless headphone has a guaranteed radio transmission throughout your entire household.

### Features

- connect to several headphones simultaneously
- radio transmission for reliable reception even through walls and ceilings of up to 80 m (freefield)
- light weight, also suitable for spectacle wearers
- transmitter, tray and charging station all in one (3-in-1)
- integrated volume- and balance controller
- optimal user and listening comfort
- automatic on/off switch
- intelligent rechargeable battery function
- nonstop operating period, approx. 4h per rechargeable battery (2nd rechargeable battery included)

### Technological Information

- wireless headphone with 3 different transmission frequency selections
- frequency response: 15–20,000 Hz
- max. sound pressure level: 120 dB
- carrier frequency: channel A 863.25 MHz, channel B 863.75 MHz, channel C 864.75 MHz
- weight (headset): approx. 59 g (incl. batteries)
- weight (transmitter): ca. 200 g
- convertible jack plug: 3.5 mm (incl. 4 adapters)

### Corporate Headquarters

audifon GmbH & Co. KG  
Werner-von-Siemens-Straße 2  
99625 Kölleda  
Germany

Phone +49-3635-4056-590  
Fax +49-3635-4056-589

[contact@audifon.com](mailto:contact@audifon.com)

[www.audifon.com](http://www.audifon.com)