

audifon
hearing systems

THAT'S HOW SMART HEARING CAN BE

●● audifon Wings hearing aids



AUDIFON WINGS - WIDE RANGE OF SMART FUNCTIONS!

We offer hearing aids which can be connected to a smartphone and to a mobile App for easy operation.

**HEARING AID FOR EVERYONE WHO
WANTS SMART CONNECTIVITY:**

- ✓ Made in Germany
- ✓ Connection to smartphone
- ✓ Make direct calls via Bluetooth
- ✓ Choose between disposable or rechargeable battery.



reddot winner 2022

LIVING
SMART CAN
BE THIS
EASY.

A woman with her back to the camera, wearing a black and white striped long-sleeved shirt, stands on a beach with her arms raised in a gesture of freedom or joy. The background shows the ocean and a clear sky.

READY
FOR THE
DAY!

POWERFUL AND COMFORTABLE THANKS TO THE LATEST BATTERY TECHNOLOGY.

Ready for the day with just 2.5 of overnight charging: lewi R Li / risa R Li are rechargeable behind-the-ear hearing aids with external receiver functionality based on our well-proven Wings chip. Enjoy the convenience of a rechargeable lithium-ion battery combined with wireless Bluetooth technology. Disposable battery is a thing of the past.

- ✓ Powerful: 28 mAh lithium-ion battery technology for up to 23 hours of wireless streaming
- ✓ Intuitive: wireless radio technology, fully integrated with Bluetooth and 2.4 GHz radio frequency
- ✓ Convenient: Pairing hearing aids with smartphone
- ✓ Sustainable: no need ever to change battery
- ✓ Versatile: Charging case doubles up as a power bank and enables up to 3 times of charging a pair of hearing aids.





MORE THAN
JUST MUSIC TO
MY EARS.

AUDIFON WINGS - MORE DIRECT LISTENING PLEASURE!

Crystal clear sound – without the annoying background noise: Switch between hearing aid mode and audio streaming mode on the easy to use audifon App.

- ✓ Answer calls directly via your audifon Wings without having to pick up the phone
- ✓ Enjoy music or audio books without the disturbing background noise
- ✓ With the **audifon multistreamer pro**, you can even stream audio from your TV directly to your hearing aid

Note: Direct streaming without accessories is available only with iOS devices and newer models of Android smartphones.
We recommend accessories for streaming from other audio sources (e.g. TV), see page 15.

TECHNOLOGY
THAT
INSPIRES

TWICE AS SIMPLE: AUDIFON WINGS HEARING AIDS.

audifon Wings is available in two attractive designs: as a powerful standard hearing aid and as a discreetly invisible hearing aid with an external receiver unit.

audifon Wings behind-the-ear hearing aid

- ✓ Simple to use and operate
- ✓ Perfect fit



audifon Wings Behind-the-ear hearing aid with external receiver unit

- ✓ Available with disposable battery or rechargeable battery technology
- ✓ Simple to use and operate
- ✓ Barely noticeable and very discreet to others
- ✓ With an even more natural sound



Both models are available in two different technology levels. Find out more on page 13.

AUDIFON WINGS - THE INTELLIGENT DAILY HELPER.

audifon Wings combine excellent sound quality with a variety of practical functions:

- ✓ **REMOTE CONTROL**
Use the audifon App as a remote control to adjust the sound, volume and much more.
- ✓ **HEARING AID LOCATION FINDER**
Easy to locate: Misplaced hearing aids can be located quickly and easily thanks to the audifon App.
- ✓ **SAVE LISTENING SITUATIONS**
Save and link frequently used settings as Favourites to a specific location. The App adds the listening programmes to Favourites and automatically activates them in the corresponding environment.
- ✓ **LIVE MICROPHONE**
Use iPhone as an additional microphone to block out distracting background noise.

Note: Android smartphones require an accessory for some functions - see page 15.

THAT’S HOW SMART
HEARING CAN BE.



AUDIFON WINGS
HEARING AIDS.

	audifon risa	audifon lewi
Rechargeable battery technology	✓	✓
audifon app	✓	✓
Audio streaming*	✓	✓
Remote control	✓	✓
Live microphone*	✓	✓
Hearing aid finder	✓	✓
Comfort features are standard for good hearing	✓	✓
Automatic recognition and adaptation to listening environment	✓	✓
Improved detection and reduction of background noise		✓
Automatic music recognition		✓
Recognition of diverse listening situations		✓
Changes to Comfort function made on one side of the hearing aid is automatically synchronised to the other side.		✓

*audifon multistreamer pro is needed for older Android smartphone models

CRYSTAL CLEAR SOUND - THANKS TO AUDIFON MULTISTREAMER PRO.

Get even more out of your audifon Wings with these perfectly matched accessories: The **audifon multistreamer pro** allows you to play external audio source directly to your hearing aids.

FOR ALL WHO APPRECIATE PERFECT SOUND:

- ✓ Crystal clear TV sound
- ✓ Stream music and audio books to your hearing aids
- ✓ Listen to conversation across a table more effortlessly
- ✓ Hands-free calling



HEARING SYSTEMS „MADE IN GERMANY“.

audifon develops and produces modern, high-performance hearing aids at its plant in Germany. Clients in over 65 countries worldwide rely on our proven hearing aid technologies.

Be convinced of the advantages of smart hearing aid technology “Made in Germany”! Get to know our Wings hearing aids.

CE 0297

The products described have undergone separate conformity assessments and fulfil the required EU regulations for medical devices. The four-digit number denotes the relevant test agency that is involved in the Conformity Assessment.

M A D E I N G E R M A N Y

