

## ●● pro care

### Easy to use cleaning products for hearing aids

With pro care, audifon offers a comprehensive range of care and cleaning products tailored to suit all audifon hearing aid types and ear moulds. Ask your hearing care professional for further information.



#### pro care cleansing wipes:

- For cleaning hearing aids and ear moulds
- Cleansing wipes moistened with cleaning solution
- Contents: 25 individually wrapped wipes



#### pro care cleaning spray:

- For quick and convenient cleaning of hearing aids and ear moulds
- With integrated cleaning brush for hard-to-reach areas on the ear mould
- Contents: 50 ml incl. spray head and brush head



#### pro care cleaning tabs:

- For deep cleansing of ear mould pieces with active oxygen
- Contents: 20 pieces individually packed



#### pro care BTE cleaning kit:

- Practical all-in-one kit for cleaning BTE devices
- Contents: 8 cleaning tabs, 1 cleaning tub and 1 drying air puffer



#### pro care cleaning tubs:

- Plastic tub with insert strainer for use with pro care cleaning tabs and water



#### Sales:

audifon GmbH & Co. KG · Werner-von-Siemens-Straße 2 · D-99625 Köllda · Germany  
Tel. +49-3635-4056-590 · Fax +49-3635-4056-589 · [contact@audifon.com](mailto:contact@audifon.com)



#### Manufacturer:

MG Développement  
[www.mg-developpement.com/](http://www.mg-developpement.com/)

## ●● pro care

### pro care disinfecting and cleansing wipes:

- Specially designed for cleaning hearing aids and ear moulds
- Cleansing wipes moistened with disinfectant cleaning solution
- Contents: Dispenser box with 25 wipes



#### Sales:

audifon GmbH & Co. KG · Werner-von-Siemens-Straße 2 · D-99625 Kölleda · Germany  
Tel. +49-3635-4056-590 · Fax +49-3635-4056-589 · [contact@audifon.com](mailto:contact@audifon.com)



#### Manufacturer:

Kleinmann GmbH  
<https://kleinmann.net/>

MG Développement  
[www.mg-development.com/](http://www.mg-development.com/)

### Reliable care products for drying hearing aids

#### pro care drying tubs:

- Plastic tubs for drying BTE and ITE hearing aids using drying capsules



#### pro care drying capsules:

- Removes condensation in BTE and ITE hearing aids to prevent damage due to moisture
- Contents: 4 drying capsules with silica gel



#### pro care drying air puffer:

- Effective removal of moisture from the hearing aid tubing and ear mould



#### Sales:

audifon GmbH & Co. KG · Werner-von-Siemens-Straße 2 · D-99625 Kölleda · Germany  
Tel. +49-3635-4056-590 · Fax +49-3635-4056-589 · [contact@audifon.com](mailto:contact@audifon.com)



#### Manufacturer:

MG Développement  
[www.mg-development.com/](http://www.mg-development.com/)