

●● sueno pro ITE



- Tinnitus-Module
- 9 WDRC-Channels
- Up to 4 Programs
- Low Battery Indicator
- Start-up Delay
- Data Logging
- nanoShield
- Volume control
- Option: Switch, Battery Size 10/312

Technical Data

Operating Voltage

EN 60118-7:2005 (2 cm³-coupler)

1.30 V

EN 60118-0/A1:1994 (Ear Simulator)

1.30 V

ANSI S3.22-2014 (2 cm³-coupler)

1.30 V

Battery Current

1.02 mA

1.02 mA

1.02 mA

Battery Type

10/312

10/312

10/312

Average Battery Life (Zinc-Air)

100/140 h

100/140 h

100/140 h

Tinnitusmasker

Noise Level (RMS)

70

79

70

Frequency Range


100 Hz–8000 Hz

200 Hz–8000 Hz

100 Hz–8000 Hz

All curves are measured with 2cm³-coupler (EN 60318-4:2010) in reference setting.

Third Octave Band Noise



Standard



Programming

- Cable: Cable set J or K
- Battery: without Battery
- Progr.-Box: HI-PRO/HI-PRO II
HI-PRO USB
NOAHlink
- Software: audifit 5.6

