

●● sueno pro S



- Tinnitus-Module
- 9 Channels
- Up to 4 Programs
- Rocker Switch (programmable)
- Low Battery Indicator
- Start-up Delay
- Data Logging
- Water repellent coating


Equipment:

- Thin Tube System

Technical Data

EN 60118-0:2015 (2 cm³-coupler) & **ANSI S3.22-2014** (2 cm³-coupler)

Operating Voltage	1,30 V
Noise Level (RMS)	69 dB
Frequency Range	100 Hz–8000 Hz
Battery Current	1,18 mA
Battery Type	312
Average Battery Life (Zinc-Air)	120 h



- 1 Rocker switch
- 2 Battery compartment/
On-Off-switch
- 3 Battery compartment lock

Standard



Programming

Cable:	Cabel Set H or I
Battery:	without Battery
Progr.-Box:	HI-PRO/HI-PRO 2 HI-PRO USB NOAHlink
Software:	audifit 5.8



●● sueno pro S

All curves are measured with 2 cm³-coupler (EN 60318-5:2015) in reference setting.

